

THERA-PAW® INC.

ASSISTIVE PRODUCTS FOR ANIMALS WITH SPECIAL NEEDS
 therapaw.com questions@therapaw.com

HINDLIMB DORSI-FLEX ASSIST

REDUCES KNUCKLING, ASSISTS W/ TARSAL FLEXION
 mobility aid for dogs with distal hindlimb weakness



Measure-ments in Inches	BOOT Measurement Range	TIBIAL CUFF Circumference Measurement A & Cuff Height	TOE STRAPS Measurement Range
	SIZE	Paw Width Paw Length tracing in standing	Circumference & Cuff Height
			From tip of toes to across back of Achilles
XX-Small <i>fits most weighing 15-30 lbs</i>	2" paw width 2 1/2" paw length <i>uses size T/B Thera-Paw Boot</i>	2 1/2" – 4" circumference Cuff height 2"	5 1/2" – 6 3/4"
X-Small <i>fits most weighing 30-40 lbs</i>	2 1/4" paw width 2 3/4" paw length <i>uses size T/C Thera-Paw Boot</i>	3 1/2" – 4 3/4" circumference Cuff height 2 1/4"	6" – 7 3/4"
Small <i>fits most weighing 35-65 lbs</i>	2 1/2" paw width 3" paw length <i>uses size T/D Thera-Paw Boot</i>	4 1/4" – 6" circumference Cuff height 3"	7 1/2" – 9"
Medium <i>fits most weighing 60-90 lbs</i>	2 3/4" paw width 3 1/2" paw length <i>uses size T/E Thera-Paw Boot</i>	5" – 6 3/4" circumference Cuff height 3 3/4"	9 1/2" – 11"
Large <i>fits most weighing 80-110 lbs</i>	3" paw width 3 3/4" paw length <i>uses size T/F Thera-Paw Boot</i>	5 1/2" – 7 3/4" circumference Cuff height 4"	10 1/2" – 13"
X-Large <i>fits most weighing 100-130 lbs</i>	3 1/2" paw width 4 1/2" paw length <i>uses size T/G Thera-Paw Boot</i>	6 1/2" – 9" circumference Cuff height 4 3/4"	11 1/2" – 14"

Notes: The Hindlimb Dorsi-Flex Assist should provide a snug fit. Tibial cuff is adjustable. Toe straps are adjustable.